



Anti-TIMM21 (aa 6-107) polyclonal antibody (DPABH-16658)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May participate in the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.
Immunogen	antigen sequence, corresponding to amino acids 6-107 of Human C18orf55.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TIMM21 translocase of inner mitochondrial membrane 22 homolog \(yeast\) \[Homo sapiens \]](#)

Official Symbol	TIMM21
Synonyms	TIMM21; translocase of inner mitochondrial membrane 21 homolog (yeast); TIM21; HSPC154; C18orf55; mitochondrial import inner membrane translocase subunit Tim21; TIM21-like protein, mitochondrial;
Entrez Gene ID	29090
Protein Refseq	NP_054896.2
UniProt ID	A8K1K8
Pathway	Metabolism of proteins;
Function	molecular_function; protein binding;