



# Anti-FAM171B (aa 491-624) polyclonal antibody (DPABH-16679)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	FAM171B (family with sequence similarity 171, member B) is a protein-coding gene. GO annotations related to this gene include DNA binding. An important paralog of this gene is FAM171A2.
<b>Immunogen</b>	antigen sequence, corresponding to amino acids 491-624 of Human FAM171B.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">FAM171B family with sequence similarity 172, member B [ Homo sapiens ]</a>
<b>Official Symbol</b>	FAM171B
<b>Synonyms</b>	FAM171B; family with sequence similarity 171, member B; KIAA1946; protein FAM171B;
<b>Entrez Gene ID</b>	<a href="#">165215</a>
<b>Protein Refseq</b>	<a href="#">NP_803237.3</a>
<b>UniProt ID</b>	<a href="#">Q6P995</a>
<b>Function</b>	DNA binding;