



Anti-NCBP2L (aa 20-122) polyclonal antibody (DPABH-16710)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Belongs to the RRM NCBP2 family. Contains 1 RRM (RNA recognition motif) domain.
Immunogen	antigen sequence, corresponding to amino acids 20-122 of Human NCBP2L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [NCBP2L nuclear cap binding protein subunit 3-like \[Homo sapiens \]](#)

Official Symbol	NCBP2L
Synonyms	NCBP2L; nuclear cap binding protein subunit 2-like; dJ820B18.1;
Entrez Gene ID	392517
Function	RNA cap binding; nucleotide binding;