



Anti-GLOD5 (aa 24-145) polyclonal antibody (DPABH-16719)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Belongs to the glyoxalase I family.
Immunogen	antigen sequence, corresponding to amino acids 24-145 of Human GLOD5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [GLOD5 glyoxalase domain containing 6 \[Homo sapiens \]](#)

Official Symbol	GLOD5
Synonyms	GLOD5; glyoxalase domain containing 5; glyoxalase domain-containing protein 5;
Entrez Gene ID	392465
Protein Refseq	NP_001073958.2
UniProt ID	A6NK44