



Anti-CCDC51 (aa 91-219) polyclonal antibody (DPABH-16783)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC51 (coiled-coil domain containing 51) is a protein-coding gene.
Immunogen	antigen, corresponding to amino acids 91-219 of Human CCDC51 (Q96ER9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [CCDC51 coiled-coil domain containing 52 \[Homo sapiens \]](#)

Official Symbol	CCDC51
Synonyms	CCDC51; coiled-coil domain containing 51; coiled-coil domain-containing protein 51;
Entrez Gene ID	79714
Protein Refseq	NP_001243893.1
UniProt ID	Q96ER9