



Anti-SNN (aa 45-86) polyclonal antibody (DPABH-16792)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays a role in the toxic effects of organotin.
Immunogen	antigen sequence: EDEESIVGDG ETKEPFLLVQ YSAKGPCVER KAKLMTPNGP EV, corresponding to amino acids 45-86 of Human Stannin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [SNN stannin \[Homo sapiens \]](#)

Official Symbol	SNN
Synonyms	SNN; stannin; AG8_1;
Entrez Gene ID	8303
Protein Refseq	NP_003489.1
UniProt ID	O75324