



Anti-LRRN1 (aa 383-508) polyclonal antibody (DPABH-16802)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LRRN1 (leucine rich repeat neuronal 1) is a protein-coding gene. Diseases associated with LRRN1 include atrioventricular septal defect , and neuronitis. An important paralog of this gene is LRFN3.
Immunogen	antigen sequence, corresponding to amino acids 383-508 of Human LRRN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRRN1 leucine rich repeat neuronal 2 [Homo sapiens]
Official Symbol	LRRN1
Synonyms	LRRN1; leucine rich repeat neuronal 1; NLRR-1; FIGLER3; leucine-rich repeat neuronal protein 1; neuronal leucine-rich repeat protein 1; fibronectin type III, immunoglobulin and leucine rich repeat domains 3;
Entrez Gene ID	57633
Protein Refseq	NP_065924.3
UniProt ID	Q6UXK5