



Anti-TMIGD2 (aa 44-146) polyclonal antibody (DPABH-16813)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMIGD2 (transmembrane and immunoglobulin domain containing 2) is a protein-coding gene. Diseases associated with TMIGD2 include bullous pemphigoid , and endotheliitis.
Immunogen	antigen, corresponding to amino acids 44-146 of Human TMIGD2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [TMIGD2 transmembrane and immunoglobulin domain containing 3 \[Homo sapiens \]](#)

Official Symbol	TMIGD2
Synonyms	TMIGD2; transmembrane and immunoglobulin domain containing 2; IGPR1; IGPR-1; transmembrane and immunoglobulin domain-containing protein 2; CD28 homolog; immunoglobulin and proline-rich receptor 1; immunoglobulin-containing and proline-rich receptor 1; transmembrane and immunoglobulin domain-containing protein 2 variant 2; transmembrane and immunoglobulin domain-containing protein 2 variant 3;
Entrez Gene ID	126259
Protein Refseq	NP_001162597.1
UniProt ID	Q96BF3
Function	coreceptor activity; protein binding;