



Anti-SMIM19 (aa 44-107) polyclonal antibody (DPABH-16835)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SMIM19 (small integral membrane protein 19) is a protein-coding gene.
Immunogen	antigen sequence: KRNKRRIMRI FSVPPTEETL SEPNFYDTIS KIRLRQQLEM YSISRKYDYQ QPQNQADSVQ LSLE, corresponding to amino acids 44-107 of Human C8orf40.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [SMIM19 small integral membrane protein 20 \[Homo sapiens \]](#)

Official Symbol	SMIM19
Synonyms	SMIM19; small integral membrane protein 19; C8orf40; UPF0697 protein C8orf40;
Entrez Gene ID	114926
Protein Refseq	NP_001129146.1
UniProt ID	Q96E16