



Anti-DYNAP (aa 143-210) polyclonal antibody (DPABH-16846)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DYNAP (dynactin associated protein) is a protein-coding gene. GO annotations related to this gene include protein binding.
Immunogen	antigen sequence: NSSIVQLST NDGECVTVKP GTPSPACPPT MTTTSTVPAS TATESTPSTA TAATTSTEPI TVAPTDHL, corresponding to C terminal amino acids 143-210 of Human C18orf26.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	DYNAP dynactin associated protein [Homo sapiens]
Official Symbol	DYNAP
Synonyms	DYNAP; dynactin associated protein; C18orf26; dynactin-associated protein; full;
Entrez Gene ID	284254
Protein Refseq	NP_775900.1
UniProt ID	Q8N1N2
Function	protein binding;