



Anti-FMR1NB (aa 98-187) polyclonal antibody (DPABH-16874)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FMR1NB (fragile X mental retardation 1 neighbor) is a protein-coding gene. Diseases associated with FMR1NB include dyslexia , and sarcoma.
Immunogen	antigen sequence: VLANGHILPN SENAHGQSLE EDSALEALLN FFFPTTCNLR ENQVAKPCNE LQDLSESECL RHKCCFSSSG TTSFKCFAPF RDVPKQMMQM, corresponding to amino acids 98-187 of Human FMR1NB.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FMR1NB fragile X mental retardation 2 neighbor [Homo sapiens]
Official Symbol	FMR1NB
Synonyms	FMR1NB; fragile X mental retardation 1 neighbor; CT37; NYSAR35; NY-SAR-35; fragile X mental retardation 1 neighbor protein; cancer/testis antigen 37; sarcoma antigen NY-SAR-35;
Entrez Gene ID	158521
Protein Refseq	NP_689791.1
UniProt ID	Q8N0W7