



Anti-HEG1 (aa 129-248) polyclonal antibody (DPABH-16911)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HEG1 (heart development protein with EGF-like domains 1) is a protein-coding gene. Diseases associated with HEG1 include cavernous malformation , and cerebral cavernous malformation. GO annotations related to this gene include calcium ion binding.
Immunogen	antigen sequence, corresponding to amino acids 129-248 of Human HEG1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	HEG1 heart development protein with EGF-like domains 2 [Homo sapiens]
Official Symbol	HEG1
Synonyms	HEG1; heart development protein with EGF-like domains 1; HEG; MST112; MSTP112; protein HEG homolog 1; HEG homolog 1; heart of glass;
Entrez Gene ID	57493
Protein Refseq	NP_065784.1
UniProt ID	Q9ULI3
Function	calcium ion binding;