



Anti-SEZ6 (aa 774-897) polyclonal antibody (DPABH-16947)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in cell-cell recognition and in neuronal membrane signaling. Seems to be important for the achievement of the necessary balance between dendrite elongation and branching during the elaboration of a complex dendritic arbor. Involved in the development of appropriate excitatory synaptic connectivity.
Immunogen	Recombinant fragment corresponding to amino acids 774-897 of Human Sez6 (UniProt: Q53EL9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SEZ6 seizure related 7 homolog (mouse) [Homo sapiens]
Official Symbol	SEZ6
Synonyms	SEZ6; seizure related 6 homolog (mouse); seizure protein 6 homolog; SEZ-6; hSEZ-6;
Entrez Gene ID	124925
Protein Refseq	NP_001092105.1
UniProt ID	Q53EL9