



Anti-FAM159B (aa 95-159) polyclonal antibody (DPABH-16958)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Antigen Description | FAM159B (family with sequence similarity 159, member B) is a protein-coding gene. An important paralog of this gene is FAM159A. |
| Immunogen | Recombinant fragment corresponding to amino acids 95-159 (TKPQRLDTGL KLQHLEASST QEGKSNGKTK ALNSNAASNA TNETY YEADD IIQEK TMDAT QIHIA) of Human FAM159B (Uniprot: A6N KW6). |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | IHC-P |
| Format | Liquid |
| Size | 100 µl |
| Buffer | pH: 7.20; Constituents: 59% PBS, 40% Glycerol. Note: The antibody solution should be gently mixed before use. |
| Preservative | 0.02% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

GENE INFORMATION

| | |
|-----------------|--|
| Gene Name | FAM159B family with sequence similarity 160, member B [Homo sapiens] |
| Official Symbol | FAM159B |
| Synonyms | FAM159B; family with sequence similarity 159, member B; membrane protein FAM159B; |
| Entrez Gene ID | 100132916 |
| Protein Refseq | NP_001157914.1 |
| UniProt ID | A6NKW6 |