



Anti-TMEM109 (aa 212-242) polyclonal antibody (DPABH-16969)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM109 (transmembrane protein 109) is a protein-coding gene. GO annotations related to this gene include molecular_function.
Immunogen	antigen sequence: LEAKVRGLER QVEELRWRQR RAAKGARSVE E, corresponding to amino acids 212-242 of Human TMEM109.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-Fr, ICC, IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM109 transmembrane protein 110 [Homo sapiens]
Official Symbol	TMEM109
Synonyms	TMEM109; transmembrane protein 109; mg23; mitsugumin-23;
Entrez Gene ID	79073
Protein Refseq	NP_076997.1
UniProt ID	Q9BVC6
Function	molecular_function;