



Anti-PRR12 (aa 966-1044) polyclonal antibody (DPABH-17016)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PRR12 (proline rich 12) is a protein-coding gene. GO annotations related to this gene include DNA binding. An important paralog of this gene is QSER1.
Immunogen	Recombinant fragment corresponding to amino acids 966-1044 (TKVKAEPKK KKKWLKEAG GNATAGGGPP GSSSDSESSP GAPSEDERAV PGRLLKTRAM REMYRSYVEM LVSTALDPD) of Human PRR12. (Uniprot:Q9ULL5)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol. Note: The antibody solution should be gently mixed before use.
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PRR12 proline rich 13 [Homo sapiens]
Official Symbol	PRR12
Synonyms	PRR12; proline rich 12; KIAA1205; proline-rich protein 12;
Entrez Gene ID	57479
Protein Refseq	NP_065770.1
UniProt ID	Q9ULL5
Function	DNA binding;