



Anti-LRRC47 (aa 163-273) polyclonal antibody (DPABH-17022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LRRC47 (leucine rich repeat containing 47) is a protein-coding gene. GO annotations related to this gene include phenylalanine-tRNA ligase activity and RNA binding. An important paralog of this gene is FARSB.
Immunogen	antigen sequence, corresponding to amino acids 163-273 of Human LRRC47.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRRC47 leucine rich repeat containing 48 [Homo sapiens]
Official Symbol	LRRC47
Synonyms	LRRC47; leucine rich repeat containing 47; leucine-rich repeat-containing protein 47;
Entrez Gene ID	57470
Protein Refseq	NP_065761.1
UniProt ID	A0A024R4G1
Function	phenylalanine-tRNA ligase activity; poly(A) RNA binding; protein binding;