



Anti-GPR22 (aa 236-308) polyclonal antibody (DPABH-17030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Orphan receptor.
Immunogen	Antigen sequence corresponding to amino acids 236-308 (QALNIRIGTR FSTGQKKKAR KKKTISLTQ HEATDMSQSS GGRNVVFGVR TSVSVIIALR RAVKRHRERR ERQ) of Human GPCR GPR22.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol. Note: The antibody solution should be gently mixed before use.
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	GPR22 G protein-coupled receptor 23 [Homo sapiens]
Official Symbol	GPR22
Synonyms	GPR22; G protein-coupled receptor 22; probable G-protein coupled receptor 22;
Entrez Gene ID	2845
Protein Refseq	NP_005286.2
UniProt ID	A4D0R8
Pathway	GPCRs, Class A Rhodopsin-like;
Function	G-protein coupled receptor activity;