



Anti-TMBIM1 (aa 25-97) polyclonal antibody (DPABH-17047)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMBIM1 (transmembrane BAX inhibitor motif containing 1) is a protein-coding gene. Diseases associated with TMBIM1 include endotheliitis. An important paralog of this gene is TMBIM4.
Immunogen	antigen: GYGQPSVLPG GYPAYPGYPQ PGYGHYPAGYP QMPPTHPMP MNYGPGHGYD GEERAVSDSF GPGEWDDRKY RHT, corresponding to amino acids 25-97 of Human TMBIM1 (Q969X1)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMBIM1 transmembrane BAX inhibitor motif containing 2 [Homo sapiens]
Official Symbol	TMBIM1
Synonyms	TMBIM1; transmembrane BAX inhibitor motif containing 1; LFG3; RECS1; MST100; PP1201; MSTP100; protein lifeguard 3; transmembrane BAX inhibitor motif-containing protein 1;
Entrez Gene ID	64114
Protein Refseq	NP_071435.2
UniProt ID	A0A024R472
Function	death receptor binding;