



Anti-C11ORF48 (aa 3-128) polyclonal antibody (DPABH-17064)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C11orf48 (chromosome 11 open reading frame 48) is a protein-coding gene.
Immunogen	antigen, corresponding to amino acids 3-128 of Human C11orf48 (Uniprot: Q9BQE6)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C11ORF48 chromosome 11 open reading frame 49 [Homo sapiens]
-----------	---

Official Symbol	C11ORF48
Synonyms	C11ORF48; chromosome 11 open reading frame 48; uncharacterized protein C11orf48;
Entrez Gene ID	79081
Protein Refseq	NP_077004.2
UniProt ID	A0A024R584