



Anti-TMEM25 (aa 262-364) polyclonal antibody (DPABH-17086)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM25 (transmembrane protein 25) is a protein-coding gene. Diseases associated with TMEM25 include breast cancer.
Immunogen	antigen, corresponding to amino acids 262-364 of Human TMEM25 (Uniprot: Q86YD3)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TMEM25 transmembrane protein 26 \[Homo sapiens \]](#)

Official Symbol	TMEM25
Synonyms	TMEM25; transmembrane protein 25; 0610039J01Rik;
Entrez Gene ID	84866
Protein Refseq	NP_001137506.1
UniProt ID	Q86YD3