



Anti-ILDR2 (aa 362-493) polyclonal antibody (DPABH-17160)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ILDR2 (immunoglobulin-like domain containing receptor 2) is a protein-coding gene. Diseases associated with ILDR2 include essential hypertension , and hypertension. An important paralog of this gene is LSR.
Immunogen	Recombinant fragment, corresponding to amino acids 362-493 of Human ILDR2 (Q71H61).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC, IB Recommended dilution: IHC: 1:50- 1:200 IB: 0.4 µg/mL
Format	Liquid
Size	100 µl
Buffer	Phosphate-buffered saline, pH 7.2, containing 40% glycerol and 0.02% sodium azide.
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ILDR2 immunoglobulin-like domain containing receptor 3 [Homo sapiens]
Official Symbol	ILDR2
Synonyms	ILDR2; immunoglobulin-like domain containing receptor 2; C1orf32; dJ782G3.1; immunoglobulin-like domain-containing receptor 2;
Entrez Gene ID	387597
Protein Refseq	NP_955383.1
UniProt ID	Q71H61