



Anti-CLMN (aa 446-553) polyclonal antibody (DPABH-17193)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CLMN (calmin (calponin-like, transmembrane)) is a protein-coding gene. Diseases associated with CLMN include emery-dreifuss muscular dystrophy , and muscular dystrophy. GO annotations related to this gene include actin binding. An important paralog of this gene is ACTN4.
Immunogen	antigen, corresponding to amino acids 446-553 of Human CLMN.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CLMN calmin (calponin-like, transmembrane) [Homo sapiens]
Official Symbol	CLMN
Synonyms	CLMN; calmin (calponin-like, transmembrane); calmin; calponin like transmembrane domain protein; calponin-like transmembrane domain protein;
Entrez Gene ID	79789
Protein Refseq	NP_079010.2
UniProt ID	Q6NUQ2
Function	actin binding;