



Rabbit anti-Human KIAA0513 polyclonal antibody (DPABH-17267)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA0513 (KIAA0513) is a protein-coding gene. Diseases associated with KIAA0513 include neuronitis.
Immunogen	antigen, corresponding to amino acids 134-271 of Human KIAA0513.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [KIAA0513 KIAA0514 \[Homo sapiens \]](#)

Official Symbol	KIAA0513
Synonyms	KIAA0513; uncharacterized protein KIAA0513;
Entrez Gene ID	9764
Protein Refseq	NP_001273494.1
UniProt ID	O60268