



Anti-IGSF10 (aa 456-601) polyclonal antibody (DPABH-17308)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in the maintenance of osteochondroprogenitor cells pool.
Immunogen	antigen, corresponding to amino acids 456-601 of Human IGSF10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [IGSF10 immunoglobulin superfamily, member 11 \[Homo sapiens \]](#)

Official Symbol	IGSF10
Synonyms	IGSF10; immunoglobulin superfamily, member 10; CMF608; immunoglobulin superfamily member 10; calvaria mechanical force protein 608;
Entrez Gene ID	285313
Protein Refseq	NP_001171616.1
UniProt ID	Q6WRI0