



Anti-GPR158 (aa 914-1052) polyclonal antibody (DPABH-17330)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Orphan receptor.
Immunogen	antigen sequence, corresponding to amino acids 914-1052 of Human GPR158.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [GPR158 G protein-coupled receptor 159 \[Homo sapiens \]](#)

Official Symbol	GPR158
Synonyms	GPR158; G protein-coupled receptor 158; probable G-protein coupled receptor 158;
Entrez Gene ID	57512
Protein Refseq	NP_065803.2
UniProt ID	Q5T848
Function	G-protein coupled receptor activity;