



Anti-PCDHB5 (aa 67-170) polyclonal antibody (DPABH-17341)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.
Immunogen	antigen sequence corresponding to amino acids 67-170 of Human PCDHB5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PCDHB5 protocadherin beta 6 \[Homo sapiens \]](#)

Official Symbol	PCDHB5
Synonyms	PCDHB5; protocadherin beta 5; PCDH-BETA5; protocadherin beta-5; PCDH-beta-5;
Entrez Gene ID	26167
Protein Refseq	NP_056484.1
Function	calcium ion binding;