



Anti-WSCD2 (aa 61-170) polyclonal antibody (DPABH-17372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	WSCD2 (WSC domain containing 2) is a protein-coding gene. Diseases associated with WSCD2 include keratosis , and hypertension. An important paralog of this gene is WSCD1.
Immunogen	antigen sequence, corresponding to amino acids 61-170 of Human WSCD2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [WSCD2 WSC domain containing 3 \[Homo sapiens \]](#)

Official Symbol	WSCD2
Synonyms	WSCD2; WSC domain containing 2; WSC domain-containing protein 2;
Entrez Gene ID	9671
Protein Refseq	NP_055468.2
UniProt ID	Q2TBF2