



Anti-RELL1 (aa 127-257) polyclonal antibody (DPABH-17393)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RELL1 (RELT-like 1) is a protein-coding gene. An important paralog of this gene is RELT2.
Immunogen	antigen, corresponding to amino acids 127-257 of Human RELT1 (UniProt ID: Q8IUW5).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [RELL1 RELT-like 2 \[Homo sapiens \]](#)

Official Symbol	RELL1
Synonyms	RELL1; RELT-like 1; RELT-like protein 1; tmp_locus_29; receptor expressed in lymphoid tissues like 1;
Entrez Gene ID	768211
Protein Refseq	NP_001078868.1
UniProt ID	Q8IUW5