



Anti-LRRC70 (aa 347-488) polyclonal antibody (DPABH-17423)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Renders cells highly sensitive to the activation by cytokines and lipopolysaccharide (LPS).
Immunogen	antigen, corresponding to amino acids 347-488 of Human LRRC70.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [LRRC70 leucine rich repeat containing 71 \[Homo sapiens \]](#)

Official Symbol	LRRC70
Synonyms	LRRC70; leucine rich repeat containing 70; SLRN; leucine-rich repeat-containing protein 70; synleurin;
Entrez Gene ID	100130733
Protein Refseq	NP_852607.3
UniProt ID	Q7Z2Q7