



Anti-SNX19 (aa 625-750) polyclonal antibody (DPABH-17454)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in several stages of intracellular trafficking.
Immunogen	antigen, corresponding to amino acids 625-750 of Human SNX19.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [SNX19 sorting nexin 20 \[Homo sapiens \]](#)

Official Symbol	SNX19
Synonyms	SNX19; sorting nexin 19; CHET8; sorting nexin-19;
Entrez Gene ID	399979
Protein Refseq	NP_055573.2
UniProt ID	B7ZAU9
Function	phosphatidylinositol binding; protein binding;