



Anti-C16ORF89 (aa 73-207) polyclonal antibody (DPABH-17465)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C16orf89 (chromosome 16 open reading frame 89) is a protein-coding gene. Diseases associated with C16orf89 include thyroiditis.
Immunogen	antigen, corresponding to amino acids 73-207 of Human C16orf89.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	C16ORF89 chromosome 16 open reading frame 90 [Homo sapiens]
Official Symbol	C16ORF89
Synonyms	C16ORF89; chromosome 16 open reading frame 89; UPF0764 protein C16orf89;
Entrez Gene ID	146556
Protein Refseq	NP_001091984.2
UniProt ID	Q6UX73
Function	protein homodimerization activity;