



Anti-NCMAP (aa 54-100) polyclonal antibody (DPABH-17484)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NCMAP (noncompact myelin associated protein) is a protein-coding gene. GO annotations related to this gene include structural constituent of myelin sheath.
Immunogen	antigen sequence, corresponding to amino acids 54-100 of Human C1orf130.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NCMAP noncompact myelin associated protein [Homo sapiens]
Official Symbol	NCMAP
Synonyms	NCMAP; noncompact myelin associated protein; MP11; C1orf130; noncompact myelin-associated protein; myelin protein of 11 kDa;
Entrez Gene ID	400746
Protein Refseq	NP_001010980.1
UniProt ID	Q5T1S8
Function	structural constituent of myelin sheath;