



Anti-LSMEM1 (aa 3-61) polyclonal antibody (DPABH-17503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LSMEM1 (leucine-rich single-pass membrane protein 1) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to N terminal amino acids 3-61 (HSSQDTGSCG IQEDGKLYVV DSINDLNKLN LCPAGSQHLF PLEDKIPVLG TNSGNGSRS) of Human C7orf53.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [LSMEM1 leucine-rich single-pass membrane protein 2 \[Homo sapiens \]](#)

Official Symbol	LSMEM1
Synonyms	LSMEM1; leucine-rich single-pass membrane protein 1; C7orf53;
Entrez Gene ID	286006
Protein Refseq	NP_001127940.1
UniProt ID	Q8N8F7