



Anti-KIAA1551 (aa 1597-1729) polyclonal antibody (DPABH-17523)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	P15RS is upregulated in cells overexpressing cyclin-dependent kinase inhibitor p15(INK4b) (CDKN2B; MIM 600431) and may have a role in cell cycle regulation (Liu et al., 2002 [PubMed 12470661]).[supplied by OMIM, Mar 2008]
Immunogen	antigen, corresponding to amino acids 1597-1729 of Human C12orf35 (Q9HCM1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA1551 KIAA1552 [Homo sapiens]
Official Symbol	KIAA1551
Synonyms	KIAA1551; C12orf35; uncharacterized protein KIAA1551; uncharacterized protein C12orf35;
Entrez Gene ID	55196
Protein Refseq	NP_060639.3
UniProt ID	Q9HCM1