



Anti-LRTM1 (aa 139-286) polyclonal antibody (DPABH-17533)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein which interacts with the N-terminus of the stem-loop binding protein (SLBP) and the 3 end of histone mRNA. This interaction facilitates the activation of histone mRNA translation. Alternative splicing results in multiple transcript variants encoding different isoforms.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 139-286 of Human LRTM1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRTM1 leucine-rich repeats and transmembrane domains 2 [Homo sapiens]
Official Symbol	LRTM1
Synonyms	LRTM1; leucine-rich repeats and transmembrane domains 1; HT017; leucine-rich repeat and transmembrane domain-containing protein 1;
Entrez Gene ID	57408
Protein Refseq	NP_065729.1
UniProt ID	Q9HBL6