



Anti-TMEM120B (aa 25-127) polyclonal antibody (DPABH-17571)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM120B (transmembrane protein 120B) is a protein-coding gene. GO annotations related to this gene include molecular_function.
Immunogen	antigen: THRIYKQKLE ELAALQTLCS SSISKQKKHL KDLKLTQRC KRHASREEAE LVQQMAANIK ERQDVFFDME AYLPKKNGLY LNLVLGNVNV TLLSNQAKFA YKD, corresponding to amino acids 25-127 of Human TMEM120B (A0PK00).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM120B transmembrane protein 121B [Homo sapiens]
Official Symbol	TMEM120B
Synonyms	TMEM120B; transmembrane protein 120B;
Entrez Gene ID	144404
Protein Refseq	NP_001074294.2
UniProt ID	A0PK00
Function	molecular_function;