



Anti-KIAA1715 (aa 254-375) polyclonal antibody (DPABH-17594)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in limb and central nervous system development.
Immunogen	antigen, corresponding to amino acids 254-375 of Human LNP (Q9C0E8).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [KIAA1715 KIAA1716 \[Homo sapiens \]](#)

Official Symbol	KIAA1715
Synonyms	KIAA1715; Ul; LNP; LNP1; ulnaless; protein lunapark; Lunapark; 2310011O18Rik; 4921514L11Rik; limb and neural patterns;
Entrez Gene ID	80856
Protein Refseq	NP_085153.1
UniProt ID	Q9C0E8