



Anti-FDCSP (aa 28-81) polyclonal antibody (DPABH-17625)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Can bind to the surface of B-lymphoma cells, but not T-lymphoma cells, consistent with a function as a secreted mediator acting upon B-cells.
Immunogen	amino acids 28-81 (KRSISDSDEL ASGFFVFPYP YPFRPLPPIP FPRFPWFRRN FPIPIPESAP TTPL) of Human FDCSP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FDCSP follicular dendritic cell secreted protein [Homo sapiens]
Official Symbol	FDCSP
Synonyms	FDCSP; follicular dendritic cell secreted protein; C4orf7; FDC-SP; follicular dendritic cell secreted peptide; FDC secreted protein;
Entrez Gene ID	260436
Protein Refseq	NP_694542.1
UniProt ID	Q540F3