



Anti-TMEM207 (aa 77-130) polyclonal antibody (DPABH-17636)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM207 (transmembrane protein 207) is a protein-coding gene. Diseases associated with TMEM207 include signet ring cell adenocarcinoma , and adenocarcinoma.
Immunogen	internal sequence amino acids 77-130 (LRRPRIDSHR RTMAVFAVGD LDSIYGTEAA VSPTVGIHLQ TQTPDLYPVP APCF) of Human TMEM207.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM207 transmembrane protein 208 [Homo sapiens]
Official Symbol	TMEM207
Synonyms	TMEM207; transmembrane protein 207; UNQ846; SRSR846;
Entrez Gene ID	131920
Protein Refseq	NP_997199.1
UniProt ID	Q6UWW9