



Anti-TMEM26 (aa 312-367) polyclonal antibody (DPABH-17646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM26 (transmembrane protein 26) is a protein-coding gene. Diseases associated with TMEM26 include alzheimers disease.
Immunogen	C terminal amino acids 312-367 (RASLRQSEG LKGEHGCRAQ TSESGPSQRD WQNESKEGLA IPLRGSPVTS DDSHHT) of Human TMEM26, isoform 1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 0.02Preservative: 0.02% Sodium azide; Constituents: 59% PBS, 40% Glycerol
Preservative	0.025% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM26 transmembrane protein 27 [Homo sapiens]
Official Symbol	TMEM26
Synonyms	TMEM26; transmembrane protein 26;
Entrez Gene ID	219623
Protein Refseq	NP_848600.2
UniProt ID	Q6ZUK4