



Anti-ZDHHC6 (aa 169-203) polyclonal antibody (DPABH-17656)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZDHHC6 (zinc finger, DHHC-type containing 6) is a protein-coding gene. GO annotations related to this gene include transferase activity, transferring acyl groups and zinc ion binding. An important paralog of this gene is ZDHHC16.
Immunogen	internal sequence amino acids 169-203 (QLYHRLSFGW NTVKIDMSAA RRDPLPIVPF GLAAF) of Human ZDHHC6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZDHHC6 zinc finger, DHHC-type containing 7 [Homo sapiens]
Official Symbol	ZDHHC6
Synonyms	ZDHHC6; zinc finger, DHHC-type containing 6; ZNF376; palmitoyltransferase ZDHHC6; DHHC-6; zinc finger protein 376; transmembrane protein H4; probable palmitoyltransferase ZDHHC6; zinc finger, DHHC domain containing 6; zinc finger DHHC domain-containing protein 6;
Entrez Gene ID	64429
Protein Refseq	NP_071939.1
UniProt ID	Q9H6R6
Function	palmitoyltransferase activity; protein-cysteine S-palmitoyltransferase activity; zinc ion binding;