



Anti-SMCO2 (aa 210-280) polyclonal antibody (DPABH-17665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SMCO2 (single-pass membrane protein with coiled-coil domains 2) is a protein-coding gene.
Immunogen	internal sequence amino acids 210-280 (LKKKVIERLE DLCKNVELLS AKLRMYQMEA EDTDSHSSEE IDTEEMEALL PQAPASFLVQ KSPPRNTAWK R) of Human C12orf70.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [SMCO2 single-pass membrane protein with coiled-coil domains 3 \[Homo sapiens \]](#)

Official Symbol	SMCO2
Synonyms	SMCO2; single-pass membrane protein with coiled-coil domains 2; C12orf70; single-pass membrane and coiled-coil domain-containing protein 2;
Entrez Gene ID	341346
Protein Refseq	NP_001138482.1
UniProt ID	A6NFE2