



Anti-OPALIN (aa 55-141) polyclonal antibody (DPABH-17674)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	OPALIN (oligodendrocytic myelin paranodal and inner loop protein) is a protein-coding gene. Diseases associated with OPALIN include fucosidosis , and temporal lobe epilepsy.
Immunogen	amino acids 55-141 (RRRSSIEAME ESDRPCEISE IDDNPKISEN PRRSPTHEKN TMGAQEAHIY VKTVAGSEEP VHDYRPTIE MERRRGLWWL VPRLSLE) of Human OPALIN.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	OPALIN oligodendrocytic myelin paranodal and inner loop protein [Homo sapiens]
Official Symbol	OPALIN
Synonyms	OPALIN; oligodendrocytic myelin paranodal and inner loop protein; TMP10; HTMP10; TMEM10; opalin; transmembrane protein 10; transmembrane protein TMP10;
Entrez Gene ID	93377
Protein Refseq	NP_001035191.1
UniProt ID	Q96PE5