



Anti-TMEM200A (aa 88-125) polyclonal antibody (DPABH-17695)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM200A (transmembrane protein 200A) is a protein-coding gene. Diseases associated with TMEM200A include medulloblastoma. An important paralog of this gene is TMEM200C.
Immunogen	fragment corresponding to internal sequence amino acids 88-125 (EHFIDAETTL STNETQVIRN EGGVVVRRFFE QHLHSDKM) of Human TMEM200A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM200A transmembrane protein 201A [Homo sapiens]
Official Symbol	TMEM200A
Synonyms	TMEM200A; transmembrane protein 200A; TTMA; TTMC; KIAA1913; two transmembrane C;
Entrez Gene ID	114801
Protein Refseq	NP_001245205.1
UniProt ID	A8K2A1