



# Rabbit Anti-Human EPGN Polyclonal Antibody (DPABH-17705)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Polyclonal to Epigen.
<b>Specificity</b>	This antibody detects endogenous levels of Epigen protein.
<b>Target</b>	Epigen
<b>Immunogen</b>	Synthesized peptide derived from human Epigen. (Internal)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µl, 200 µl
<b>Buffer</b>	PBS containing 50% glycerol and 0.5% BSA
<b>Preservative</b>	0.02% Sodium Azide

**Storage**

Store at -20°C, and avoid repeat freeze-thaw cycles.

---

## BACKGROUND

**Introduction**

Promotes the growth of epithelial cells. May stimulate the phosphorylation of EGFR and mitogen-activated protein kinases.

---

**Keywords**

EPGN;epithelial mitogen;EPG;PRO9904;ALGV3072;epigen;epithelial mitogen homolog;

---

## GENE INFORMATION

**Entrez Gene ID**

[255324](#)

---

**UniProt ID**

[Q6UW88](#)

---