



Anti-PCNXL2 (aa 341-447) polyclonal antibody (DPABH-17716)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in tumorigenesis of colorectal carcinomas with high microsatellite instability (MSI-H).
Specificity	DPABH-17716 shows specific reactivity against isoform 1 of PCNXL2.
Immunogen	internal sequence amino acids 341-447 of Human PCNXL2, isoform 1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PCNXL2 pecanex-like 3 (Drosophila) [Homo sapiens]
Official Symbol	PCNXL2
Synonyms	PCNXL2; pecanex-like 2 (Drosophila); pecanex-like protein 2;
Entrez Gene ID	80003
Protein Refseq	NP_055616.3
UniProt ID	A6NKB5