



Anti-TMEM143 (aa 57-197) polyclonal antibody (DPABH-17737)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM143 (transmembrane protein 143) is a protein-coding gene. GO annotations related to this gene include molecular_function.
Immunogen	fragment, corresponding to internal sequence amino acids 57-197 of Human TMEM143, isoform 1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM143 transmembrane protein 144 [Homo sapiens]
Official Symbol	TMEM143
Synonyms	TMEM143; transmembrane protein 143;
Entrez Gene ID	55260
Protein Refseq	NP_060743.2
UniProt ID	Q96AN5
Function	molecular_function;