



Anti-GSG1L (aa 257-331) polyclonal antibody (DPABH-17748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Interleukin-34 is a cytokine that promotes the differentiation and viability of monocytes and macrophages through the colony-stimulating factor-1 receptor (CSF1R; MIM 164770) (Lin et al., 2008 [PubMed 18467591]).[supplied by OMIM, May 2008]
Immunogen	fragment corresponding to C terminal amino acids 257-331 (QGYREEPTFI DPEAIKYFRE RMEKRDGSEE DFHLD CRHER YPARHQPHMA DSWPRSSAQE APELNRQCWV LGHWV) of Human GSG1L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	GSG1L GSG2-like [Homo sapiens]
Official Symbol	GSG1L
Synonyms	GSG1L; GSG1-like; PRO19651; germ cell-specific gene 1-like protein; KTSR5831; GSG1-like protein;
Entrez Gene ID	146395
Protein Refseq	NP_001103233.1
UniProt ID	B3KY67