



Anti-TMEM223 (aa 119-202) polyclonal antibody (DPABH-17770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM223 (transmembrane protein 223) is a protein-coding gene.
Immunogen	fragment corresponding to amino acids 119-202 (RSVRSVVLRA GGQQVTLTTH APFGLGAHFT VPLKQVSCMA HRGEVPAMLP LKVKGRRFYF LLDKTGHFPN TKLFDNTVGA YRSL) of Human TMEM223.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM223 transmembrane protein 224 [Homo sapiens]
Official Symbol	TMEM223
Synonyms	TMEM223; transmembrane protein 223;
Entrez Gene ID	79064
Protein Refseq	NP_001073970.1
UniProt ID	A0PJW6